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NOKANO | **USER MANUAL**

TOOL DISINFECTION TUB

BEFORE USE, READ THE USER MANUAL AND KEEP THIS USER MANUAL

PRODUCT PURPOSE

The tray is intended for:
1. Disinfection of metal and stainless-steel tools,
2. Use exclusively with disinfectant liquids intended for cosmetic tools.

The manufacturer is not responsible for any damage resulting from improper use or failure to follow the instructions provided in the manual.

PRODUCT SPECIFICATION

1. Material: Plastic
2. Dimensions: 225 x 120 x 60 mm

INSTRUCTIONS FOR USE

1. Select a disinfectant suitable for the intended purpose.
2. Prepare the solution according to the disinfectant manufacturer's instructions.
3. Open the lid of the tray.
4. Pour the prepared disinfectant solution into the tray.
5. Clean tools of any biological residue before disinfection.
6. Place them in the inner basket of the tray.
7. Fully immerse the tools in the solution.
8. Follow the exposure time specified by the disinfectant manufacturer.
9. Remove the basket with the tools.
10. Rinse the tools with clean water and dry thoroughly.
11. Dispose of the used solution in accordance with the disinfectant manufacturer's instructions.
12. After each use, wash the tray and dry it thoroughly.

PRECAUTIONS

1. This product is not a toy – keep out of reach of children.
2. Use only on stable, flat surfaces.
3. Do not use aggressive chemicals that may damage the plastic.
4. Avoid skin contact with disinfectant solutions – wear protective gloves.
5. Do not expose the tray to high temperatures (› 50°C).
6. Wash the product with warm water and detergent before first use.
7. Do not use in a dishwasher or microwave oven.

STORAGE CONDITIONS

1. Store in a dry, shaded place.t
2. Protect from direct sunlight and heat sources.

DISPOSAL

1. Dispose of in accordance with local plastic waste regulations.
2. The product does not require special treatment as hazardous waste.

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